

Deliver Quality Physio Treatment in Home Healthcare

# Physio Kit Us FDA Cleared Devices for Quickly and Conveniently deliver the best treatment in pair and in the convenient of the pair and the convenient of th

promote tissue healing at HOME



# Stim<sup>3</sup>

Prove your professional practice with most efficient technology in hand and deliver IFT, TENS & EMS.



4 channels



pole application



5000 Hz Frequency



LCD Digital display



Operated



100% modulation

# JUS-1

Easy to use 1 MHz ultrasound provides continuous and pulsed treatment; Highly effective to increase local blood flow. reduce swelling, chronic inflammation and promote healing process.





Microprocessor Controlled Unit




Continuous and Pulsed Mode



Digital Display of Selections



... deliver the quality treatment at home



## **Gate Control Protocol / Opiate Protocol**

Immediate and Lasting Analgesic effect in chronic & acute pain conditions, to facilitate functional movement during ADLs.

#### **Trigger Point Protocol**

Deactivate TP activity to facilitate functional movement

# **CHF / Endurance Protocol**

Increase endurance of muscle group to improve ADLs, stability, postural tone.

## **Explosive Protocols**

Increase contraction power and speed of muscle groups to improve stability.

#### **Functional Exercise**

Increase ability of muscle group to contract and improve muscle recruitment.

## **Sports Medicine**

Musculoskeletol problems like inflammation from injuries, muscles sprains or runners knee. Increase the stretchiness of muscles and healing rate.

Pre-operative and post-operative physiotherapeutic intervention



- 1800 102 8684 (Toll Free)
- mktg@joharidigital.com
- www.joharidigital.com